



HiSense™

Mycoplasma PCR Detection Kit

Version:6.0

Latest update: 2024.02.01

Instruction Manual

Cat. No. HD-100

Research Use Only. Not for Use in Diagnostic Procedures.

Contents

1. Product Information	3
Introduction	
Kit Specificity	
Kit Sensitivity	
Prevention of carryover contamination	
Precaution for cross contamination	
2. Contents and Storage	4
Materials Provided	
Storage	
3. Test protocol	5
Workflow	
Prepare the Template (Sample)	
Prepare for PCR	
4. Results / Trouble Shooting	8
Gel Electrophoresis Result	

1. Product Information

Introduction

HiSense™ Mycoplasma PCR Detection Kit utilizes the PCR, which is established as the method of choice for highest sensitivity in the detection of mycoplasma, ureaplasma and acholeplasma contamination in cell culture and other cell culture derived biologicals. The primers are specific to the highly conserved 16S ribosomal RNA sequences in the mycoplasma genome. This allows for detection of all mycoplasma species including *Acholeplasma laidlawii*, *M. arginini*, *M. fermentans*, *M. gallisepticum*, *M. genitalium*, *M. hominis*, *M. hyorhinae*, *M. orale*, *M. pneumoniae*, *Spiroplasma citri* and *U. urealyticum*. All mycoplasma species including 8 genus and about 209 species can be detected simultaneously. The “European Pharmacopoeia” and Guideline of MFDS recommend checking for unspecific detection of *Clostridium*, *Lactobacillus*, and *Streptococcus*. The 16S ribosomal RNA of other bacteria such as *E. coli*, *Clostridium*, *Lactobacillus*, *Streptococcus*, and plant and animal cells are not amplified.

Kit Specificity

The HiSense™ Mycoplasma PCR Detection Kit detects Mycoplasma species simply, reliably, and rapidly. To detect the presence of these microorganisms, the assay uses the polymerase chain reaction (PCR) to amplify a target unique to a wide variety of mycoplasmas. The kit can detect more than 209 different Mycoplasma species, including *Acholeplasma laidlawii* and *Spiroplasma citri*. The kit does not detect other genera or cell-line DNA.

Kit Sensitivity

The sensitivity of the PCR using this kit is 1 to 10 copies of the target DNA per reaction. Sensitivity of the assay in real culture samples depends on the quality of the sample preparation.

Prevention of carryover contamination

Take precautions to minimize the potential for carryover of nucleic acids from one experiment to the next. Use separate work areas and pipettes for template preparation and PCR setup steps.

The 2X PCR Master Mix contains dUTP instead of dTTP. When dUTP replaces dTTP in PCR amplification, UNG treatment (Uracil-N-glycosylase, HiSense™ UDG heat labile Cat.No.SCU-100, not provided in this kit) can prevent the subsequent reamplification of dU-containing PCR products. UNG acts on single- and double-stranded dU-containing DNA by hydrolysis of uracil-glycosidic bonds at dU-containing DNA sites. With this strategy, carryover contamination is eliminated, while the template DNA (DNA containing T) remains intact.

Precaution for cross contamination

Cross contamination may occur while handling the positive control DNA contained in the kit. To prevent cross contamination, use separate locations and pipettes for handling positive control during the PCR process. Be sure to use PCR grade water included in the kit and wear gloves during experiments. Additionally, use filter tips to prevent aerosols.

2. Contents and Storage

Materials Provided

Label	Cap	HD-100
2X PCR Master Mix	Blue	1.5 ml
Primer Mix	Red	200 µl
Positive Control DNA	Yellow	50 µl
DNase Free Water	White	1.2 ml

Storage

Upon receipt, store at -20°C.

Note

- 1) Repeat thawing reduces quality of product.
- 2) If frequent freeze and thaw is needed, aliquot the products and use in order.
- 3) Check the label on the product for details.

3. Test protocol

Workflow


Sample Preparation

A gDNA Extraction Method



Mycoplasma harvest

- 1) Transfer 1 ml of cell culture supernatant.
- 2) Centrifuge (12,000~15,000 rpm, 10min) and discard supernatant.



Genomic DNA extraction

- 3) Follow the DNA extraction kit manual.
- 4) 50µl of the elution volume is obtained in the elution step

B Boiling Method



Remove cell debris

- 1) Transfer 1 ml of cell culture supernatant.
- 2) Centrifuge (1,000~1,500 rpm, 5min) and transfer 1ml of supernatant to a new tube.



Mycoplasma harvest & washing

- 3) Centrifuge (12,000~15,000 rpm, 10min) and discard supernatant.
- 4) Add 1 ml D.W or Tris-EDTA buffer and repeat step 3.



Boiling

- 5) Discard supernatant, add D.W 100 µl and vortexing.
- 6) Boiling (98°C, 10 min)
- 7) Centrifuge (12,000~15,000 rpm, 5 min)
- 7) Transfer 90 µl of supernatant to a new tube.

Reagent Preparation



2X PCR Master Mix



Primer Mix



Positive Control DNA




DNase Free Water



Template DNA

	Sample	PC	NTC
2X PCR Master Mix	15µl	15µl	15µl
Primer Mix	2µl	2µl	2µl
Positive Control DNA	-	1µl	-
Template DNA	5µl	-	-
DNase Free Water		up to 30µl	

PCR Setup and Run



Analyze the results



Prepare the Template (Sample)

Samples should be derived from cultures which are at 90-100% confluence. Penicillin and streptomycin in the culture media do not inhibit PCR or affect test sensitivity. To avoid false positive results, we recommend the use of the PCR grade water delivered with the kit, aerosol-preventive filter tips and gloves.

[A] Genomic DNA Extraction Method

- 1) Transfer 1 ml of cells with media to a 1.5 ml clean tube.
- 2) Centrifuge the tube at 12,000~15,000 rpm for 10 minutes to pellet mycoplasma and discard supernatant.
- 3) From the next step, follow the DNA extraction kit manual you are using.
- 4) In the elution step, elution is performed with 50 µl, and the resulting genomic DNA is used as a PCR test sample.

[B] Boiling Method

- 1) Transfer 1 ml of cell culture supernatant to a 1.5 ml clean tube.
- 2) Centrifuge at 1,000~1,500 rpm for 5 minutes to pellet cell debris and transfer 1 ml of supernatant to a clean tube.
- 3) Centrifuge the tube at 12,000~15,000 rpm for 10 minutes to pellet mycoplasma and discard supernatant. Mycoplasma pellet may not be seen in naked eyes.
- 4) Add 1 ml DNase free water or Tris-EDTA buffer for washing and repeat step 3).
This process must be performed so that the PCR reaction is not inhibited.
- 5) Discard supernatant completely, add DNase free water 100 µl into mycoplasma pellet and mix with vortexing.
- 6) Heat the supernatant at 98°C for 10 minutes and centrifuge at 12,000~15,000 rpm for 5 minutes.
Caution!! Be careful when you heat the sample at 98°C. Heating it in PCR machine with heating cover.
- 7) Transfer 90 µl of supernatant to a clean tube. 5 µl of supernatant is used as the template in PCR.

Prepare for PCR

1. Prepare the set of reactions listed in the following table.

(Caution!! Don't vigorous vortexing.)

Reaction Components		Sample Reaction	Positive Control Reaction	No Template Control Reaction
Label	Cap			
2X PCR Master Mix	Blue	15 µl	15 µl	15 µl
Primer Mix	Red	2 µl	2 µl	2 µl
Template DNA	-	5 µl	-	-
Positive Control DNA	Yellow	-	1 µl	-
DNase Free Water	White	Up to 30 µl		
Final volume		30 µl	30 µl	30 µl


2. Set up the qPCR instrument to run the PCR cycling (amplification) program specified below.

Steps & Cycle		Temp(°C)	Time
Pre-denature		95	5 min
PCR	40 Cycles	Denature	95 20 sec
		Anneal	60 30 sec
		Extend	72 30 sec

4. Results / Trouble Shooting

Gel Electrophoresis Result

Tested with Quantitative *M. arginini* Genomic DNA (Cat. No. QAGD, CellSafe)



1. 10^4 copies
2. 10^3 copies
3. 10^2 copies
4. 10^1 copies
5. 1 copy
6. Positive Control
- 7,8. No Template Control
9. 100bp ladder

When mycoplasma contamination exists, a band with around 250-300bp appears. An internal control DNA band with around 700bp means the right performance of PCR reaction.

Using the following table, determine whether the test cell culture is infected with mycoplasma.

Mycoplasma Band	Internal Control DNA Band	Interpretation
Positive	Positive	Mycoplasma contamination
Positive	Negative*	Mycoplasma contamination
Negative	Positive	Mycoplasma non-contamination
Negative	Negative	PCR inhibition or inadequate sample preparation

* In case of severe mycoplasma contamination, Internal Control DNA band can be not detected.

- 1) Recommend performing no template control reaction and positive control reaction by adding 1 μ l of positive control DNA.
- 2) If the PCR reaction is inhibited by high FBS concentration, the use of genomic DNA as a template may be helpful.
- 3) PCR inhibiting substances may accumulate in the medium of cells (e.g., Hybridoma cell, etc.). In this case, the use of diluted sample or genomic DNA as a template may be helpful.